

Microplastic Polymer Analysis Report

Report ID	e8166b93-a2a6-4b13-930e-f892bc23de7e
Generated On	20 May 2026 11:39
Status	Completed
Sample ID	HDPE004

Prediction Results

Detected Polymer	HDPE
Confidence Score	100.00%
Model Used	Peak Identifier

FTIR Spectral Analysis

- 2915 cm^{-1} → CH_2 asymmetric stretching
- 1470 cm^{-1} → CH_2 bending

Confidence Interpretation

■■■■■■■■■■ 100.00%

Scientific References

[1] Harper, C. Handbook of Plastics Technologies.

[2] Thompson, R. Plastic debris in oceans.

[3] Barnes, D. Accumulation of plastics.

[4] Andrady, A. Environmental impacts of plastics.

[5] Jung, M. FTIR analysis of HDPE.

AI Generated Summary

The analyzed sample was identified as High-Density Polyethylene (HDPE) with a confidence score of 100.00%. The FTIR spectral characteristics strongly support the presence of HDPE polymer chains. The detected CH_2 asymmetric stretching at 2915 cm^{-1} and CH_2 bending at 1470 cm^{-1} are characteristic of the hydrocarbon functional groups in HDPE. This polymer is widely used in plastic bags, bottles, and containers. The analysis indicates a reliable classification suitable for environmental assessment and waste management studies.