

Microplastic Polymer Analysis Report

Report ID	dc18df42-71e8-4067-8f9b-bee6ab862547
Generated On	18 May 2026 18:42
Status	Completed
Sample ID	HDPEFTIR001

Prediction Results

Detected Polymer	HDPE
Confidence Score	100.00%
Model Used	Peak Identifier

FTIR Spectral Analysis

- 2915 cm^{-1} $\rightarrow$ CH₂ asymmetric stretching
- 1470 cm^{-1} $\rightarrow$ CH₂ bending

Confidence Interpretation

■■■■■■■■■■ 100.00%

Scientific References

- [1] Harper, C. Handbook of Plastics Technologies.
- [2] Thompson, R. Plastic debris in oceans.
- [3] Barnes, D. Accumulation of plastics.
- [4] Andrady, A. Environmental impacts of plastics.
- [5] Jung, M. FTIR analysis of HDPE.

AI Generated Summary

The analyzed sample was identified as HDPE with a confidence score of 100.00%. The FTIR spectral characteristics strongly support the presence of HDPE polymer chains. The detected functional groups and

absorption bands are consistent with known spectral signatures of HDPE. The analysis indicates a reliable classification suitable for environmental and material identification studies.